

2~8°C /-10~-26°C Frigorífico y congelador combinados	
Modelo	YCD-EL450
Tipo de gabinete	Vertical
Capacidad (L)	450, R:225, F:225
Tamaño interno (An. x Pr. x Al.) mm	R:650*570*627, F:650*570*627
Tamaño externo (An. x Pr. x Al.) mm	810*735*1960
Tamaño del paquete (An. x Pr. x Al.) mm	895*820*2135
NW/GW(Kgs)	146/156
Performance	
Temperature Range	R:2~8, F:-10~-26
Ambient Temperature	16-32°C
Cooling Performance	R:5, F:-26
Climate Class	N
Controller	Microprocessor
Display	Digital display
Refrigeration	
Compressor	2pcs
Cooling Method	R: Forced air cooling, F:Direct cooling
Defrost Mode	R:Automatic, F:Manual
Refrigerant	R600a
Insulation Thickness(mm)	R:80, F:80
Construction	
External Material	PCM
Inner Material	Aluminum plate with spraying
Shelves	R:3(coated steel wired shelf), F:6(ABS)
Door Lock with Key	Yes
Lighting	LED
Access Port	2pcs. Ø 25 mm
Casters	4 (2 casters with brake)
Data Logging/Interval/Recording Time	USB/Record every 10 minutes / 2 years

Backup Battery	Yes
Door with Heater	Yes
Alarm	
Temperature	High/Low temperature, High ambient temperature
Electrical	Power failure, Low battery
System	Sensor failure, Door ajar, Built-in USB datalog failure, Communication failure
Accessories	
Standard	RS485, contacto de alarma remota

